

ramé-hart instrument co.

19 Route 10 East, Suite 11
Succasunna, NJ 07876 USA
USA





telephone [973-448-0305](tel:973-448-0305)
fax line [973-448-0315](tel:973-448-0315)
email carl@ramehart.com
web www.customspinnerets.com





ramé-hart Prebuilt Customizable Linear Spinneret Array Kit p/n 100-10-PCLSAK

What's in the box

- (1) Stainless steel array body
- (7) Male Luer fittings for needles (output)
- (1) Female Luer fitting for syringe (inlet)
- (11) Plug screws
- (19) Sealing washers
- (1) Allen key (for plug screws)

Many configurations (shown below) are possible with 10, 15, 20, 30, 40, 50 mm output spacing. Flexible inlet mounting.





For more information on this and other spinneret products, visit:

www.customspinnerets.com

ramé-hart instrument co.